ISOLITE® ISODRY® ISOVAC®

Mouthpieces

Instructions for Use	EN
Brugsanvisning	DA
Gebrauchsanweisung	DE
Οδηγίες χρήσης	EL
Instrucciones de uso	ES
Käyttöohjeet	FI
Mode d'emploi	FR
Istruzioni per l'uso	IT
Gebruiksaanwijzing	NL
Bruksanvisning	NO
Instrukcja użtykowania	PL
Instruções de utilização	PT
Bruksanvisning	SV
Kullanım Talimatları	TR

Instructions for Use: Isolite® Mouthpieces	EN
Description	
Non-sterile Mouthpieces developed by Zyris, Inc. are intended for use with the Isolite®, Isodry® and Isovac®. Mouthpieces are single use only and available in various sizes. Mouthpieces are constructed of a soft, flexible polymer and are not made with natural rubber latex. They do not contain phthalates or BPAs. For more information on sizes, go to Zyris.com.	
How Supplied	
Mouthpieces (sold separately) are supplied in convenient operatory boxes. Each Mouthpiece is individually packaged. Non-sterile Mouthpieces are packaged inside a polyethylene bag.	
Contraindications	
No known contraindications. Do not use without reviewing how to select the correct Mouthpiece size for each patient. Do not use if unable to position the Mouthpiece correctly.	
Side Effects	
Use of the Mouthpiece may trigger a pharyngeal reflex (laryngeal spasm) in hypersensitive patients, or if it is not placed correctly within the oral cavity.	
Selecting the Mouthpiece	
To select the correct Mouthpiece size, measure the patient's interincisal opening. If between sizes, you may select the smaller size. For additional assistance, watch a video demonstration online at Zyris.com.	
Approximate Interincisal Opening of Mouthpiece Sizes	
Pediatric (PED): < 30 mm	
Extra Small (XSM): < 30 mm	
Small (SML): 30 mm	
Medium (MED): 30-45 mm	
Medium Deep Vestibule (MDV): 30-45 mm	
Large (LRG): > 45 mm	
Attaching the Non-Sterile Mouthpiece	
Slide the Mouthpiece onto the device until the two locking barbs engage. Water or mouthpiece lubricant may be used to lubricate the inside surface of the Mouthpiece sleeve to ease attachment.	
Placing the Mouthpiece	
Prior to placing the Mouthpiece within the patient's oral cavity, moisten the Mouthpiece and the patient's lips.	
Step 1	
Hold the sides of the Control Head with your index finger and thumb. With your other hand, fold the Cheek Shield onto the Tongue Retractor. Gently slide the folded mouthpiece into the buccal vestibule on the side to receive treatment. Angle the mouthpiece to allow the lower edge of the Tongue Retractor to move along the buccal edge of the teeth.	
Step 2	
Move the bite block onto occlusal surfaces of the teeth, just distal to the mandibular cuspid. Instruct patient to "rest gently" on the bite block to secure.	
Step 3	
Place the Cheek Shield into the buccal vestibule. Move the Isthmus onto the retromolar pad behind the maxillary tuberosity. Adjust the Tongue Retractor in the lingual vestibule as needed. Move the bite block distally to provide more vertical working room.	
WARNINGS	
WARNING: Mouthpieces are single use only. Reuse can result in the transmission of infectious pathogens. Mouthpieces must be discarded after use. Mouthpieces cannot be safely sterilized after use.	
WARNING: Do not use if the individual Mouthpiece package has been previously damaged or opened. Non-sterile Mouthpieces are individually packaged in polyethylene bags.	
WARNING: Be sure to use established protocols and patient monitoring when using the Mouthpiece in conjunction with sedation dentistry. Patient management by a certified anesthesiologist is strongly recommended.	
PRECAUTIONS	
CAUTION: When performing dental procedures requiring high levels of blood evacuation, it is recommended that the Mouthpiece be intermittently flushed with water to prevent coagulation in the Mouthpiece to avoid loss of suction and light transmission.	

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ZYRIS™

Emergo Europe
Prinsesegracht 20
2514 AP The Hague
The Netherlands

Australian Sponsor
Emergo Australia
Level 20
Tower II, Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

96-0473 Rev F

Zyris, Inc.
6868A Cortona Drive
Santa Barbara, CA 93117 USA
Phone +1-(805) 560-9888
Fax +1-(805) 966-6416
For the USA, toll free:
1-800-560-6066
www.Zyris.com

© 2018 Zyris, Inc.
All rights reserved.
Patent Information
www.Zyris.com/patents

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