Instructions for Use

EN

Brugsanvisning

DA

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DE

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ISOLITE® ISODRY® ISOVAC®

Mouthpieces

Better Isolation = Better Dentistry®

Instructions for Use: Isolite® Mouthpieces

Description

Non-sterile and sterile Mouthpieces developed by Isolite Systems are intended for use with the Isolite®, Isodry®, and Isovac®. Mouthpieces are single use only and available in various sizes. Mouthpieces are constructed of a soft, flexible polymer and are not made with natural rubber latex. They do not contain phthalates or BPA's. For more information on sizes, go to isolitesystems.com.

How Supplied

Mouthpieces (sold separately) are supplied in convenient operatory boxes. Each Mouthpiece is individually packaged. Non-sterile Mouthpieces are packaged inside a polyethylene bag while Sterile Mouthpieces are packaged inside a thermoform tray with a Tyvek lid. Sterile Mouthpieces are sterilized using e-beam radiation. When sterilized, yellow indicator dots on both the Tyvek lid and the box turn red. Do not use the Sterile Mouthpiece if the indicator dots are not red.

Contraindications

No known contraindications. Do not use without reviewing how to select the correct Mouthpiece size for each patient. Do not use if unable to position the Mouthpiece correctly.

Side Effects

Use of the Mouthpiece may trigger a pharyngeal reflex (laryngeal spasm) in hypersensitive patients, or if it is not placed correctly within the oral cavity.

Selecting the Mouthpiece

To select the correct Mouthpiece size, measure the patient's interincisal opening. If between sizes, you may select the smaller size. For additional assistance, watch a video demonstration online at isolitesystems.com.

Approximate Interincisal Opening of Mouthpiece Sizes

Pediatric (PED): < 30 mm

Extra Small (XSM): < 30 mm

Small (SML): 30 mm

Medium (MED): 30-45 mm

Medium Deep Vestibule (MDV): 30-45 mm

Large (LRG): > 45 mm

Attaching the Non-Sterile Mouthpiece

Slide the Mouthpiece onto the device until the two locking bars engage. Water or mouthpiece lubricant may be used to lubricate the inside surface of the Mouthpiece sleeve to ease attachment.

Attaching the Sterile Mouthpiece

Follow the established protocol in your practice for handling sterile product. Open the thermoform tray by removing the Tyvek lid and place the Mouthpiece in a sterile field. Slide Mouthpiece onto the device until the two locking bars engage. Sterile water may be used to lubricate the inside surface of the Mouthpiece sleeve to ease attachment.

Placing the Mouthpiece

Prior to placing the Mouthpiece within the patient's oral cavity, moisten the Mouthpiece and the patient's lips.

Step 1

Hold the sides of the Control Head with your index finger and thumb. With your other hand, fold the Cheek Shield onto the Tongue Retractor. Gently slide the folded mouthpiece into the buccal vestibule on the side to receive treatment. Angle the mouthpiece to allow the lower edge of the Tongue Retractor to move along the buccal edge of the teeth.


Step 2

Move the bite block onto occlusal surfaces of the teeth, just distal to the mandibular cuspid. Instruct patient to "rest gently" on the bite block to secure.

Step 3

Place the Cheek Shield into the buccal vestibule. Move the Isthmus onto the retromolar pad behind the maxillary tuberosity. Adjust the Tongue Retractor in the lingual vestibule as needed. Move the bite block distally to provide more vertical working room.

WARNINGS

⚠️ **WARNING:** Mouthpieces are single use only. Reuse can result in the transmission of infectious pathogens. Mouthpieces must be discarded after use. Mouthpieces cannot be safely resterilized after use.

⚠️ **WARNING:** Do not use if the individual Mouthpiece package has been previously damaged or opened. Non-sterile Mouthpieces are individually packaged in polyethylene bags. Sterile Mouthpieces are individually packaged in a thermoform tray with a sealed Tyvek lid.

⚠️ **WARNING:** Be sure to use established protocols and patient monitoring when using the Mouthpiece in conjunction with sedation dentistry. Patient management by a certified anesthesiologist is strongly recommended.

PRECAUTIONS

⚠️ **CAUTION:** Do not use the Sterile Mouthpiece if the sterilization indicator dot on the Tyvek lid is not red.

⚠️ **CAUTION:** When performing dental procedures requiring high levels of blood evacuation, it is recommended that the Mouthpiece be intermittently flushed with water to prevent coagulation in the Mouthpiece to avoid loss of suction and light transmission.

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All rights reserved.
Innerlite, Inc.