

14 Million Pediatric Patients Protected

ZYRIS®

The **all-in-one Isolite Isolation System** allows you to achieve optimal visibility, moisture and aerosol control, patient compliance, improved efficiency and consistent results – while providing superior patient comfort and safety!



"I could not do what I do as a pediatric dentist without it."

- Dr. Anna Fencel

*Isolite has a unique ability to **settle younger patients** when placed:*

- No taste of fluids or nasty etch
- Continuous evacuation, minimizing fear of drowning or swallowing fluid
- Child has control of fluid with their tongue by moving it to the base of the bite block where it is evacuated.
- Won't slip out once placed




Isolite Pro with Posterior Mouthpiece

*An Isolite System delivers **advanced isolation and access** with:*

- True Shadowless Illumination for YOU and your TEAM – See More, Miss Less!
- Two Quadrant Isolation for faster and efficient procedures
- Dual-Channel Suction Control with continuous hands-free saliva and water evacuation
- Tongue and cheek retraction for ideal access and control
- Integrated Bite block and airway protection for patient comfort and safety

You'll love our Pediatric Mouthpiece Sizes - a Perfect Fit Every Time!



6 size options for a truly comfortable patient experience with every procedure

Studies have shown that an Isolite System can dramatically cut procedure time and boost the productivity of the entire office, allowing you to perform additional procedures and increase profitability.

- ✓ Sealants
- ✓ Restorations
- ✓ Stainless Steel Crowns
- ✓ Extractions
- ✓ Operating Room Procedures (featuring our new Operating Room Adapter)
- ✓ Pulpotomies/Pulp Therapies

Maximize your Time in the Operating Room *with the Zyris Isolite Operating Room Adapter*

The **NEW** biodegradable Operating Room Adapter allows the use of the Isolite Posterior Mouthpiece in an Operating Room setting. The single-use adapter enables quick operating room turnover time by allowing everything from the vacuum canister to the mouthpiece to be disposed of.

"I would not practice in the operating room without it. It makes practicing sedation dentistry faster, easier, and safer. It suctions all the fluids, I don't need to use a rubber dam, and I can go work on an entire side at a time."

- Dr. Ana Vázquez



Discover why over 2,600 pediatric dentists use Zyris Isolite dental isolation in their practice.

To place an order, visit zyris.com,
email us at orders@zyris.com
or call us at 1.800.560.6066

To schedule a **FREE DEMO**
visit ZyrisNow.com or
call Zyris at 1.800.560.6066